



**Many possibilities
during cataract surgery.
More possibilities after.**

Knowing your options can help
you focus on the breathtaking
experiences still to come.

So you've got cataracts. What does that mean?

For one, it means you're not alone. Currently, **about 20.5 million Americans have a cataract** in at least one eye. Cataracts are **a natural part of the aging process**.

Cataracts develop gradually, and eventually they can make everyday life difficult. It can seem as if you're viewing the world through a foggy window. You may experience:

- **Blurred vision**
- **Glare**
- **Sensitivity to light**
- **Difficulty reading**

Cataracts: A fact of life

- Usually start to develop around age 40
- No. 1 cause of impaired vision in the U.S.
- Cannot be prevented

Simulated image
with cataracts

Simulated image
without cataracts



You have an opportunity to change how you see the world

Imagine being able to see more than just blurry shapes in the distance. Good news! With cataract surgery, you can.

About **4 million procedures** take place each year in the U.S. to help patients rediscover the beauty of a clear horizon.

Cataract surgery is one of the safest and most effective surgeries in the U.S.: 98% of all cataract surgeries are performed each year without serious complications being reported.

What's more, cataract surgery gives you the potential to **treat cataracts, astigmatism and presbyopia** at the same time (if desired), **which may reduce the need for glasses.**

What's your best option?

Keep reading for insight into:

- Cataracts, astigmatism and presbyopia
- Your treatment options
- Working with your eye care specialist to choose the right option



Understanding your eyes

As with cataracts, there's nothing you can do to prevent astigmatism and presbyopia. However, you may be able to **treat them at the same time you treat your cataracts**. The treatment you choose could leave you with **better vision than you had before your cataracts**!


A closer look at eyesight

Vision is all about light rays entering and passing through the eye:


1. When you look at an object, light rays reflect off that object and reach the **cornea** on the surface of your eye.
2. The **cornea** bends the light and sends it to the **lens** inside your eye.
3. The lens further bends the light, sending it to the **retina**.
4. The **retina** interprets the light as an image and sends it to your brain.



Healthy Eye



Eye with Cataract



What is a cataract?

As you age, the lens becomes cloudy, making it harder for light to reach the retina. **This clouding of the natural lens is what we call a cataract.**



Healthy lens



Cataract lens

What is astigmatism?

In astigmatic eyes, the cornea is unevenly curved, distorting light rays before they reach the lens.

Astigmatism is present from birth in many cases. **If you have significant astigmatism, you've likely worn glasses most of your life.**



*Normal corneas are spherical, like a basketball.
Astigmatic corneas are shaped more like a football.*

What is presbyopia?

Even if you've never heard the word "presbyopia," you're probably familiar with the condition.

Presbyopia is the reason you need reading glasses when you pick up a newspaper or a restaurant menu.

Presbyopia occurs when the lens in your eye stiffens with age, losing the ability to focus light to see things close up.

So how does cataract surgery

Your surgeon will **remove the cloudy lens** in your eye and replace it with a new lens (known as an IOL). There have never been more choices than now.

Choosing the IOL that Best Fits Your Lifestyle

AcrySof® IQ ReSTOR®+2.5 D Multifocal Toric IOL with ACTIVEFOCUS™ Optical Design

Astigmatism-Correcting Multifocal Lens

- **Distance vision:** Clear
- **Astigmatism correction:** Yes
 - You may enjoy less need for prescription glasses or contacts
- **Presbyopia correction:** Yes
 - Reduced dependence on reading glasses
 - Offers quality vision for daily activities like computer work and driving

AcrySof® IQ ReSTOR®+2.5 D IOL with ACTIVEFOCUS™ Optical Design

Multifocal Lens

- **Distance vision:** Clear
- **Astigmatism correction:** No
 - Prescription glasses or contacts still needed
- **Presbyopia correction:** Yes
 - Reduced dependence on reading glasses
 - Offers quality vision for daily activities like computer work and driving

AcrySof® IQ Toric IOL

Astigmatism-Correcting Lens

- **Distance vision:** Clear
- **Astigmatism correction:** Yes
 - You may enjoy less need for prescription glasses or contacts
- **Presbyopia correction:** No
 - Reading glasses still needed

AcrySof® IQ IOL

Monofocal Lens

- Most basic cataract surgery option
- **Distance vision:** Clear
- **Astigmatism correction:** No
 - Prescription glasses or contacts still needed
- **Presbyopia correction:** No
 - Reading glasses still needed

work?

eye, then **replace it with an artificial lens**






choices when it comes to the IOL you receive.

Near
vision

Intermediate
vision

Distance
vision

Astigmatism



For more information about cataracts and cataract surgery options:

- Visit www.mycataracts.com
- Call 1-844-MY-CATARACT
- Contact your eye care specialist

ACRYSOF® IQ RESTOR® FAMILY OF MULTIFOCAL IOLS IMPORTANT PRODUCT INFORMATION FOR PATIENTS

CAUTION: Restricted by law to sale by or on the order of a physician. **DESCRIPTION:** The AcrySof® IQ ReSTOR® multifocal Intraocular Lenses (IOLs) are artificial lenses implanted in the eye of adult patients following cataract surgery. These multifocal lenses (sometimes called “presbyopia correcting IOLs”) are designed to allow for clear distance, intermediate, and near vision with the potential to be more independent of the need to use glasses for daily tasks. The AcrySof® IQ ReSTOR® Family includes the AcrySof® IQ ReSTOR® multifocal Toric IOLs which correct for astigmatism as well as cataracts. **WARNINGS / PRECAUTIONS:** You may experience and need to contact your eye doctor immediately if you have any of the following symptoms after cataract surgery: itching, redness, watering of your eye, sensitivity to light. The safety and effectiveness of AcrySof® IQ ReSTOR® IOLs have not been established in patients with eye conditions, such as an increase in eye pressure (glaucoma) or complications of diabetes in the eye (diabetic retinopathy). As with any surgical procedure, there are risks involved. These risks may include but are not limited to infection, damage to the lining of the cornea, the retinal layer which lines the inside back wall of your eye may become separated from the tissue next to it (retinal detachment), inflammation or swelling inside or outside the eye, damage to the iris (the colored diaphragm around the pupil), or an increase in eye pressure. With the AcrySof® IQ ReSTOR® multifocal IOLs, there may be a loss of sharpness of your vision that may become worse in dim light or in foggy conditions. There is also a possibility that you may have some visual effects such as rings or circles (halos) around lights at night. You may also have trouble seeing street signs due to bright lights or glare from oncoming headlights. **ATTENTION:** As with any surgical procedure, there are risks involved. Prior to surgery, ask your eye doctor to provide you with the Patient Information Brochure for the proposed AcrySof® IQ ReSTOR® multifocal IOL to be implanted. This document will further inform you of the risks and benefits associated with this IOL. Discuss any questions about possible risks and benefits with your eye doctor, your medical condition and any eye disease you may have.

ACRYSOF® IQ INTRAOCULAR LENS IMPORTANT PRODUCT INFORMATION CAUTION: Restricted by law to sale by or on the order of a physician. **DESCRIPTION:** The AcrySof® IQ Intraocular Lenses (IOLs) are artificial lenses implanted in the eye of adult patients following cataract surgery. These lenses are designed to allow for clear distance vision. However, you will likely still need glasses for reading and for distance vision particularly if you already have astigmatism. **WARNINGS/PRECAUTIONS:** You may experience and need to contact your eye doctor immediately if you have any of the following symptoms while using the antibiotic eye drops prescribed by your doctor: itching, redness, watering of your eye, sensitivity to light. The safety and effectiveness of the AcrySof® IQ IOL has not been established in patients with certain eye conditions, such as an increase in eye pressure (glaucoma) or complications of diabetes in the eye (diabetic retinopathy). As with any surgical procedure, there are risks involved. These risks may include, but are not limited to, infection, damage to the ocular structures: lining (inner surface) of the cornea, damage to the iris (the colored diaphragm around the pupil), the retinal layer which lines the inside back wall of your eye may become separated from the tissue next to it (retinal detachment), inflammation or swelling inside or outside the eye, an increase in eye pressure that may not be controlled by medicine and secondary surgical procedure. There is a possibility that this IOL could be placed incorrectly or could move within the eye. This may result in less improvement or a reduction in vision, or it may cause visual symptoms. **ATTENTION:** As with any surgical procedure, there are risks involved. Prior to surgery, discuss any questions you may have about possible risks and benefits with your eye doctor.

ACRYSOF® IQ TORIC IOL IMPORTANT PRODUCT INFORMATION CAUTION: Restricted by law to sale by or on the order of a physician. **DESCRIPTION:** The AcrySof® IQ Toric Intraocular Lenses (IOLs) are artificial lenses implanted in the eye of adult patients following cataract surgery. These lenses are designed to correct pre-existing corneal astigmatism, which is the inability of the eye to focus clearly at any distance because of difference curvatures on the cornea, and provide distance vision. **WARNINGS/PRECAUTIONS:** Contact your eye doctor immediately if you have any of the following symptoms while using the antibiotic eye drops prescribed by your doctor: itching, redness, watering of your eye, sensitivity to light. The safety and effectiveness of the AcrySof® IQ Toric IOL has not been established in patients with eye conditions, such as an increase in eye pressure (glaucoma) or complications of diabetes in the eye (diabetic retinopathy). As with any surgical procedure, there are risks involved. These risks may include, but are not limited to, infection, damage to the lining of the cornea, the retinal layer which lines the inside back wall of your eye may become separated from the tissue next to it (retinal detachment), inflammation or swelling inside or outside the eye, damage to the iris (the colored diaphragm around the pupil), an increase in eye pressure that cannot be controlled by medicine and secondary surgical procedure. A toric IOL corrects astigmatism only when it is placed in the correct position in the eye. There is a possibility that the toric IOL could be placed incorrectly or could move within the eye. This may result in less improvement or a reduction in vision because your astigmatism has not been fully corrected, or it may cause visual symptoms. **ATTENTION:** As with any surgical procedure, there are risks involved. Prior to surgery, ask your eye doctor to provide you with an AcrySof® IQ Toric Patient Information Brochure, which will inform you of the risks and benefits associated with this IOL. Discuss any questions about possible risks and benefits with your eye doctor.

